



### Interfering Substances












The following substances, naturally present in respiratory specimens or that may be artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. None of them were found to affect test performance of The COVID-19 Antigen Rapid Test Device.

Substance	Concentration	Substance	Concentration
3 OTC nasal sprays	10%	Guaiacol glycerol ether	20 mg/mL
3 OTC mouthwashes	10%	Mucin	1%
3 OTC throat drops	10%	Mupirocin	250 µg/mL
4-acetamidophenol	10 mg/mL	Oxymetazoline	25 µg/ml
Acetylsalicylic acid	10 mg/mL	Phenylephrine	10 mg/mL
Albuterol	10 mg/mL	Phenylpropanolamine	1 mg/mL
Chlorpheniramine	5 mg/mL	Zanamivir	10 mg/mL
Dexamethasone	50 µg/ml	Adamantanamine	500 ng/mL
Dextromethorphan	10 µg/ml	Oseltamivir phosphate	10 mg/mL
Diphenhydramine	5 mg/mL	Tobramycin	10 mg/mL
Doxylamine	1 mg/mL	Triamcinolone	14 mg/mL
Flunisolide	25 µg/ml		

### LITERATURE REFERENCES

1. Forni, D., Cagliani, R., Clerici, M. & Sironi, M. Molecular evolution of human coronavirus genomes. Trends Microbiol. 25, 35–48 (2017).
2. Ithete, N. L. et al. Close relative of human Middle East respiratory syndrome coronavirus in bat, South Africa. Emerg. Infect. Dis. 19, 1697–1699 (2013).

### GLOSSARY OF SYMBOLS

	Catalog number		Temperature limitation
	Consult instructions for use		Batch code
	In vitro diagnostic medical device		Use by
	Manufacturer		Contains sufficient for <n> tests
	Do not reuse		Authorized representative in the European Community
	CE marking according to IVD Medical Devices Directive 98/79/EC		



Assure Tech. (Hangzhou) Co., Ltd.  
Building 4, No. 1418-50, Moganshan Road,  
Gongshu District, Hangzhou,  
310011 Zhejiang, P.R.China

Lotus NL B.V.  
Koningin Julianaplein 10, le Verd,  
2595AA, The Hague, Netherlands

